

Bubble CPAP Circuit



CT-19250

It comprises of the following items -

- WT-96622 Air-Bubbler
- WT-96622 Heated Wire Circuit for Inflow
- HC-73138 Disposable Humidification Chamber, Auto Feeding
- CT-13540 Pressure Manifold

- Compatible to F&P MR-850 Humidifier
- Spiral Heated Wire Circuit to ensure minimum condensation during oxygen therapy (NCPAP)
- Nasal Tubing - 70 mm included in the pack.

Air-Bubbler provides a proven method of delivering 2-10 cm H₂O with gas flow from 1-12 Lpm.




Air Bubbler



Specification

Gas Inlet Connector	15M/22F connector
Gas Flow Range	4-15 LPM
Adjustable Range	2-10 cm H ₂ O
Accuracy	±1 cm H ₂ O
Water Capacity	400 ml

CT-19240



Pressure Manifold

Relief Valve with T-Piece connector & tube to connect with air blender
12 LPM maximum air inlet